



Dead Gene Walking

Team Prasad-Feldgarden:

Arjun Prasad, Team co-Lead

Michael Feldgarden, Team co-Lead

Adrien Assie, Tech Lead

Ana Ramos, Writer



EB Dickinson

Chienchi Lo

Erin Young

Breaking a gene can lead to resistance

- **Genes that are the target of antibiotics**
 - *K. pneumoniae* CirA - siderophore receptor used to import iron and the antibiotic cefiderocol
- **Porins that let antibiotics into the cell**
 - *K. pneumoniae* OmpK35/K36
 - Loss-of-function mutations confer resistance to carbapenems and other beta-lactams
- **Regulatory genes**
 - *P. aeruginosa* NalD - efflux system regulator
 - *P. aeruginosa* AmpD - beta-lactamase gene regulator
 - confer resistance to carbapenems and other beta-lactams
 - *A. baumannii* AdeS - efflux pump regulator




Project Overview


- Existing AMR gene detection software screens for known mutations
- Develop a tool and database to identify loss-of-function mutations in genes known to confer resistance when disrupted
 - Identify lesions that should 'break' the gene:
 - Frameshifts
 - Internal stops (nonsense mutations)
- Run the tool(s) on assemblies in NDARO to characterize the presence of novel loss-of-function mutations in those genes for three taxa
 - *Acinetobacter baumannii*
 - *Klebsiella pneumoniae*
 - *Pseudomonas aeruginosa*

What are Nonsense Mutations and Frameshifts?

NON SENSE MUTATION



FRAMESHIFT



Dead Gene Scanner.R




- Use a curated database of reference genes
- Blastn reference gene in whole genomes
 - Aim to find genes that missed annotation prediction
 - Find stop mutation in existing genes
- Checks and reports for STOP point mutation or frameshift in results
 - Check 1000 genomes/min



Technical approaches

- **Assembly sequence to identify broken genes**
- Translated alignment with BLASTX
 - Stop codon detection
 - A python script to run BLASTX and check the ***** position in the alignment.
 - Frameshift detection
 - “Diamond” blastx has option **--frameshift** to allow frameshift gap aligned and show **/** or **** in the alignment result
 - In-frame mutation is noted with **-** in the alignment. However, it is common and may not cause genes broken. (default detection is disabled.)

```
usage: DGW_blast.py [-h] -i INPUT_FASTA -d DATABASE [--cov COV] [--id ID] [--indels] [--diamond]
                    [--version]
```



```
Detect nonsensus/frameshift mutations by running blastx on contig fasta against target protein database
```



options:

-h, --help	show this help message and exit
-i INPUT_FASTA, --input_fasta INPUT_FASTA	contig fasta file
-d DATABASE, --database DATABASE	protein (diamond) blast DB
-o OUTPUT, --output OUTPUT	output tsv file
--cov COV	target coverage
--id ID	hit identity
--indels	check in-frame insertion and deletion
--diamond	use diamond blastx
--threads THREADS	threads
--verbose	Show more information in log
--debug	Show more information in log
--version	show program's version number and exit

Test Analyses - AST testset





- 928 isolates
 - 589 with resistant phenotype
 - 493 mutations found
 - 130 mutations were present in MicroBIGG-E too
 - 363 new mutations in 353 isolates
 - 339 with sensitive phenotype



Test Analyses - Random testset

- 3000 isolates for *Acinetobacter baumannii*, *Pseudomonas aeruginosa* and *Klebsiella pneumoniae*
 - DGW identified 877 mutations
 - 150 mutations were present in MicroBIGG-E too
 - 727 new mutations in 353 isolates



Conclusions and Future Directions

- We identified a novel OmpK35 frameshift mutation in isolate CP021955.1/SAMN04014948 which is highly carbapenem resistant and has no carbapenemases
- Compare results from the two tools on a large set of assemblies
- Run the tools on entire NDARO dataset to characterize frequency of these mutations
 - 106,000 *Klebsiella pneumoniae*
 - 39,000 *Pseudomonas aeruginosa*
 - 36,000 *Acinetobacter baumannii*